



DOCTOR MANUAL

A comprehensive guided explanation of all doctor functionalities and system processes necessary to make proper use of this unit and make a proper evaluation of patients based on available information.

CREDITS

MOTIVSMA | Carmen Areses, Pilar Bourg, Pablo Hervalejo, Carlota Laverón, Alejandra O'Shea

Table of contents:

1. Introduction	3
1.1. Role.....	3
2. Objectives	4
2.1. Primary	4
2.2. Secondary	4
3. Getting started	5
3.1. How to initialize the application.....	5
3.1.1. Through Docker Deployment.....	5
3.1.2. Manual setup	5
3.2. The database	6
3.2.1. Connecting the database	7
3.3. Register.....	8
3.4. Log in.....	9
5. Patient handling and data visualization	11
5.1. Accept/reject patient requests	11
5.2. See patients	12
5.3. See measurement sessions	12
5.4. Generate Report (PDF).....	13
5.5. See all generated reports.....	14
6. Logout.....	15

1. Introduction

The Doctor application provides medical professionals with a interface to manage their assigned patients, review pending patient requests, access medical information, and interact with the measurement sessions stored in a database.

1.1. Role

The role of the Doctor application is to enable healthcare professionals to:

- Review their assigned patients
- Inspect patient information (demographics, sessions, measurements)
- Approve or reject patient–doctor association requests
- Access physiological measurements and symptoms retrieved from the database.

All this information will allow the doctor to assess the evolution of Spinal Muscular Atrophy (SMA) patients at a distance through this telemedicine application.

2. Objectives

2.1. Primary

The main goal of this unit is to provide doctors with an interface so they can access measurement sessions sent by patients to maintain a periodic control over the patient's evolution during treatment stages.

Doctors may access their patients and sent measurement sessions. These include and electrocardiogram (ECG), electromyogram (EMG) and symptoms.

2.2. Secondary

It was considered important to provide an intuitive and easy-to-use interface which allows them to travel throughout the unit.

Additional actions (see *Section 4*) were added to ensure organization in doctor-patient assignation and to facilitate evolution following for the doctors.

3. Getting started

You can start the platform in two different ways. The first option is to use Docker, which provides a straightforward setup with all services pre-configured and ready to run. This approach is recommended if you want quick and consistent deployment across different environments without the need to have pre-requisites installed.

Alternatively, you can download the individual repositories and run each component separately. This method offers more flexibility for development and debugging, as it allows you to modify or replace specific modules without affecting the entire platform.

3.1. How to initialize the application

3.1.1. Through Docker Deployment

Docker Deployment provides a fully containerized environment for running the SMA (Spinal Muscular Atrophy) telemedicine system. In this section how to use and install it is explained below. Nevertheless, for more information refer to the README section in the GitHub repository.

1. Make sure Docker has been installed in your computer through the official website: <https://www.docker.com>
2. Download the GitHub repository on your computer terminal:

<https://github.com/pilarbourg/telemedicine-deploy>:

```
git clone https://github.com/pilarbourg/telemedicine-deploy
```

From this point on, please make sure you have all necessary requirements specified in the README and the certificate in place (Refer to the Certificate Manual)

3. Decompress the zip file and open a terminal or console on your device and navigate to the project directory.
4. Start all services by running the following in your terminal:

```
docker-compose up -d
```

5. The web app should now be accessible locally at <https://127.0.0.1>

Finally, to shut down all containers run the following in the same terminal: `docker-compose down`


3.1.2. Manual setup

In order to access the webpage:

1. Download the GitHub repository: <https://github.com/alejandraoshea/sma-server>
2. Download the GitHub repository: <https://github.com/alejandraoshea/sma-client>
3. Please ensure you have certificate (*refer to the Certificate Manual*)
4. Locate the index.html (`/sma-client/index.html`) and double click on it to open the browser

idea	30/11/2025 11:44	Carpeta de archivos	
.vscode	23/11/2025 13:13	Carpeta de archivos	
frontend	29/11/2025 22:33	Carpeta de archivos	
out	20/11/2025 16:14	Carpeta de archivos	
src	29/11/2025 13:58	Carpeta de archivos	
.DS_Store	28/11/2025 18:04	Archivo DS_STORE	11 KB
.gitignore	23/11/2025 13:13	txtfile	1 KB
index.html	23/11/2025 13:13	Chrome HTML Docu...	8 KB
README.md	5 11:04	Archivo MD	3 KB
sma-patient.iml	5 22:26	Archivo IML	1 KB

Now, you should see the initial register and log-in dashboard:




3.2. The database

The architecture of the database is designed around three main roles and their interactions: Patient, Doctor and Users and the measurement data generated by the patients:

- All authentication (Log-in/Sign-in) information is stored in the `app_users` table, which defines login credentials and role (ADMIN, DOCTOR, PATIENT) for each account.
- The patients and doctors extend the user records by storing their corresponding user information in and linking back to `app_users` through a shared `user_id`. Patients can be assigned to doctors and doctors may supervise multiple patients (for more information refer to the DOCTOR MANUAL and/or USER MANUAL).
- To support the data recollection by the patients, a `measurement_session` table which represents individual measurement events performed by a patient at a specific timestamp, and signals, which store the processed physiological data recorded during those sessions (EMG or ECG).

- Each signal entry references both the patient and the session it belongs to, allowing full traceability of all collected data.



3.2.1. Connecting the database

To connect the database stored in your computer to the server you will need to have the DUMP file that was given with the documentation, it must follow exact template on *section 5.1*.

Now that the database is created and stored in your computer you will have to go into *src/main/resources/application-local.yml* and change the following parameters:

```
spring:
  config:
    activate:
      on-profile: local

  datasource:
    url: jdbc:postgresql://localhost:5432/database_name
    driver-class-name: org.postgresql.Driver
    hikari:
      schema: public

server:
  port: 8443
  ssl:
    enabled: true
    key-store: /path/to/your/keystore.p12
    key-store-password: YOUR_KEYSTORE_PASSWORD
    key-store-type: PKCS12
    key-alias: YOUR_KEY_ALIAS

admin:
```

username: ADMIN_USERNAME
password: ADMIN_PASSWORD

operator:
username: OPERATOR_USERNAME
password: OPERATOR_PASSWORD

jwt:
secret: JWT_SECRET
expiration: 3600000

Highlights must correspond to your admin parameters; these will be used to log into the system

3.3. Register

If you don't have an account, please complete the registration form with the following required information:

1. Email: Enter your email address (this will be used as your identifier).
2. Password: Enter a password (please remember the password for future use).
3. Role: Select Doctor from the dropdown menu. You must select one option.

MotivSMA
Our website is dedicated to providing information, resources, and support for individuals and families affected by SMA. We aim to raise awareness about the condition, promote early diagnosis, and offer guidance on managing the disease.

Through our telemedicine platform, users can upload their current biosignals via a BitAlino sensor, connect with healthcare professionals, and analyze their results.

Services

Information Hub
Comprehensive resources about SMA, including symptoms, diagnosis, treatment options, and care strategies.

Telemedicine Platform
Upload biosignals using BitAlino sensors and connect with healthcare professionals for remote consultations and analysis.

Community Support
Forums and support groups for individuals and families to share experiences, advice, and encouragement.

Get Started
Whether you're a patient or doctor, sign up today to get started with MotivSMA.

Register

Email:

Password:

Select Role:

Doctor

MotivSMA
About Privacy Policy Contact
© 2023 MotivSMA. All rights reserved.

Finally, click the registration button.

MotivSMA
Our website is dedicated to providing information, resources, and support for individuals and families affected by SMA. We aim to raise awareness about the condition, promote early diagnosis, and offer guidance on managing the disease.

Through our telemedicine platform, users can upload their current biosignals via a BitAlino sensor, connect with healthcare professionals, and analyze their results.

Services

Information Hub
Comprehensive resources about SMA, including symptoms, diagnosis, treatment options, and care strategies.

Telemedicine Platform
Upload biosignals using BitAlino sensors and connect with healthcare professionals for remote consultations and analysis.

Community Support
Forums and support groups for individuals and families to share experiences, advice, and encouragement.

Get Started
Whether you're a patient or doctor, sign up today to get started with MotivSMA.

Register

Email:

Password:

Select Role:

Register


MotivSMA
About Privacy Policy Contact
© 2023 MotivSMA. All rights reserved.

3.4. Log in


Please insert your personal information after the first Log-in.

You must log-in into the doctor unit using the *email* and the *password* specified during the “Sign-in” following these steps:

- Email: Enter the email address you used during registration.
- Password: Enter your password.
- Action: Click the "*Login*" button.




Whenever these parameters are correct, you will have access to the doctor dashboard.



4. Profile Configuration

Once logged in, you can edit your profile information.

- Access the main screen (*Doctor Dashboard*).
- Locate the "Edit profile" section in the top left corner.



The image shows two UI components. On the left is the 'Edit Profile' form, which includes input fields for 'Name' and 'Surname', a 'Gender' dropdown menu with the text 'Select gender', and a 'Locality' dropdown menu with the text 'Select locality'. A red 'Save' button is at the bottom. On the right is a calendar grid for November 2025, with dates from 28/11/2025 to 28/11/2025 (repeated for each day of the month).

- The following personal details can be changed:
 - a) **Name**
 - b) **Surname**
 - c) **Gender**
 - d) **Location**

Click on “Save” to update the profile in the system.

5. Patient handling and data visualization

5.1. Accept/reject patient requests

Once you are logged into the doctor unit, the first thing you will see is the section named “*Patient requests*”. In this section, a list of pending requests from patients is available; if no patients have sent requests, this section will be empty and titled “*Patient requests (0)*”.

The screenshot shows the Doctor Dashboard interface. On the left, there is a 'Current Profile' section with a red circular icon containing a white heart and a person silhouette. Below the icon, the profile information is displayed: Name: Julian, Surname: Alvarez, Gender: MALE. On the right, there is a 'Patient Requests (0)' section, which is currently empty. Below it, there is a 'Current Patients' section containing a table with patient information and a list of sessions.

ID	Name	Surname	Gender	Birth Date	Height (cm)	Weight (kg)	Sessions
10	Diego	Maradona	MALE	1950-11-28	170	85	View Sessions
9	Mario	Kempes	MALE	1954-11-01	185	75	Hide Sessions

Below the table, there is a list of sessions:

- ▶ 28/11/2025, 18:02:15
- ▶ 28/11/2025, 16:05:17
- ▶ 28/11/2025, 14:38:48
- ▶ 28/11/2025, 13:49:34

Now you will be able to accept or decline these requests. If a patient is declined it will disappear from this section, but do not worry, they will be able to send another request to you or another doctor. If the patient request is accepted, all patient information will be displayed in a list containing all your patients.


The screenshot shows the Doctor Dashboard interface after a patient request has been accepted. The 'Patient Requests (1)' section now displays a request from Franco Colapinto, with details: Gender: MALE, Height: 175 cm, Birthdate: 2003-05-27, Weight: 68 kg. There are 'Approve' and 'Reject' buttons, with 'Approve' circled in red. Below this, the 'Current Patients' table is updated to include the new patient. The 'Patient Reports' section is also visible at the bottom.

ID	Name	Surname	Gender	Birth Date	Height (cm)	Weight (kg)	Sessions
10	Diego	Maradona	MALE	1950-11-28	170	85	View Sessions
11	Kylian	Mbappe	MALE	1998-12-20	178	75	View Sessions
9	Mario	Kempes	MALE	1954-11-01	185	75	View Sessions

5.2. See patients


As a doctor you will see a list of all your patients in a section titled “*Your Patients*”. For the list to be revealed you will have to click the “*See all patients*” button in the doctor dashboard. The patient list includes the following information of each person:

- ID
- Name
- Surname
- Gender
- Birth date
- Height
- Weight
- A button to access measurement sessions sent by the patient



5.3. See measurement sessions

Once a patient has sent their recorded signals and symptoms, you will be able to access all this information by clicking the button “*View sessions*” at the right of the information row for each patient.



Once this button is clicked a list of dropdowns with all measurement sessions registered by that specific patient will be visible. Measurement session will be named after the time stamp they were taken at.

Once a measurement session dropdown is clicked, all registered symptoms and both signals (ECG and EMG) for that time stamp can be seen.


The screenshot shows a user interface for patient management. On the left, there is a patient profile form with fields for Name (Julian), Surname (Alvarez), Gender (MALE), and Locality (Madrid). Below this is an 'Edit Profile' button. To the right, a measurement session is selected for the date 28/11/2025 at 16:05:17. This session details are enclosed in a red bracket. The details include a list of symptoms: MUSCLE_ATROPHY, MUSCLE_WEAKNESS, and TREMORS. Below the symptoms are two signal graphs: an EMG Signal and an ECG Signal. At the bottom of the session details is a text area for 'Doctor's notes here ...' and a blue 'Create Report' button.


5.4. Generate Report (PDF)

This feature allows doctors to export a formal clinical document for a specific session, useful for medical records. The same reports will be visible for your clients.

1. Navigate to a Session: Follow the steps to "See Measurement Sessions" for a patient.
2. Select Session: Click on a specific date (► 12/11/2024, 24:09:58) to expand the details.
3. Generate: Locate and click the button labelled "*Create Report*".
 - You can add any

This screenshot is identical to the one above, showing the patient profile and measurement session details. The 'Create Report' button at the bottom of the session details is highlighted with a red oval, indicating the next step in the process.

- *System Action:* The server compiles the patient's data, the list of symptoms, the signal graphs and comments you can add to the report into a single PDF file.
4. Download: The browser will automatically download a file named report_[reportId].pdf.
- Now the PDF file will appear in your downloads in the top-right section of your browser: 




- We recommend saving the report in your local computer and changing the name, so it is easier to identify where the report came from. For example: *MarioKempes9_12/11/2024Moring.pdf*
- Now your patient will be able to see their own generated reports with your notes.

Note: Please consider that if sensitive information can be inferred from the report and must be written in the notes, we recommend contacting the patient before saving and uploading the report into the server.

5.5. See all generated reports

As a doctor you are able to see all report generated by you in the bottom section of your dashboard called “Patient Reports”



6. Logout

To ensure the security of medical data, it is mandatory to close the session when you finish using the platform.

1. The navigation bar at the top of the screen.
2. Click the "Logout" button in the top right corner.

